










RoboICSI[®] Processor Icons and their functions

<i>Icons</i>	<i>Functions</i>
Mount 	To change/remove the RoboICSI [®] Holder
View 	Gets the RoboICSI [®] Holder in the same plane as the specimen 
Grasp 	Allows the RoboICSI [®] Joystick to hold the specimen in place <i>Move joystick Right to grasp</i>  <i>Move joystick Left to Release</i>
OFF 	Closes an ongoing function and enables the transition onto another.
Capture Image 	Takes a snapshot of the display screen
DaQ 	Logs Pre/During/Post-ICSI details
Shutdown 	Shuts down the RoboICSI [®]